


Biochemical characterisation


Biochemical Laboratory/Quality Control - Bioteck S.p.A.

Quantity of residual DNA



The DNA content is directly proportional to the cell population present in the tissue. The "control" is untreated pericardium. The other columns show the constant reduction of residual DNA following the first treatment cycle according to the Zymo-Teck® process. Further cycles help guarantee complete deantigenation.

Quantification of the lipid content



The "control" is untreated pericardium. The other columns show the drastic reduction of the lipid content after the first treatment cycle according to the Zymo-Teck® process.

Preserved collagen

MWM + Heart®



Electrophoretic separation on denaturing gel (SDS-PAGE). The first column (MWM) shows the molecular weight markers, the second (+), a type I collagen sample treated with pepsin, an enzyme able to cut the proteins at specific points. The third column shows a sample Heart® membrane subjected to the same treatment. The bands visible in the column (+) are entirely comparable with those seen in the Heart® membrane, confirming the presence of collagen in its native conformation.

BiOTECK®

The science of bone tissue

Bioteck S.p.A.

Headquarters:

Via E. Fermi 49 - 36057 Arcugnano (Vicenza) - Italy
Tel. +39 0444 289366 - fax: +39 0444 285272
info@bioteck.com - www.bioteck.com

Production and R&D Center:

Via G. Agnelli, 3 - 10020 Riva presso Chieri (Turin) - Italy



Headquarters

Bioteck® is an Italian company producing bone substitutes and protective membranes that are successfully used in orthopaedics, neurosurgery, oral and maxillofacial surgery.

Founded in 1995, the company continues to grow constantly and now operates in more than 50 countries around the world. A firm commitment to scientific research forms the basis for the innovative solutions offered by Bioteck® products.



Production and R&D Center

The company collaborates on numerous national and international research projects, which have driven the basic research and helped in writing important chapters in bone biology.

The in-depth knowledge acquired by Bioteck® through its research ensures the absolute quality of its products, which are subjected to strict environmental and quality controls, thereby guaranteeing a product meeting the highest quality and safety standards. Bioteck® applies a policy of total transparency, opening up the doors of its Production and R&D Center for the monitoring of its innovative manufacturing process and the intense scientific research carried out by its staff.


Quality and safety guarantee



Biochemical Laboratory/Quality Control

Heart® - an equine pericardium membrane range for neurosurgery, orthopaedics, general surgery, oral and maxillofacial surgery.


Heart® and Zymo-Teck® are registered Bioteck S.p.A. trademarks



concept: mauro forlani - vi - rev. 150.09_0_HR_C_EN01-00

BiOTECK®

The science of bone tissue



completely biocompatible

hard-wearing and suitable for suture

naturally re-absorbable

Heart® Patch

decellularized biological matrix



Heart® Patch

decellularized biological matrix

Zymo-Teck® Process: the secret of the quality of grafts

Heart® Patch is a biological matrix for dural replacement and repair, obtained from carefully-selected equine pericardium, to ensure just the right quality and thickness, and which is thereafter treated using the exclusive **Zymo-Teck® process**.

This sophisticated Bioteck® patented biochemical processes enables the elimination of all potentially immunogenic elements: cells, lipids and noncollagenous proteins, combining oxidative phases with the use of enzymes, natural catalytic proteins. The low temperatures at which the Zymo-Teck® process operates and a careful selection of the enzymes used, enables the **collagen fibres** to be preserved in their **natural conformation**.

Sterilisation with beta rays guarantees the complete elimination of any pathogens in respect of the three-dimensional structure of the fibres, keeping the excellent product **biomechanical characteristics** unaltered.

Heart® Patch acts as an **inert, biological scaffold**, able to support the tissue repair process without causing undesired reactions in the surrounding tissues. Once applied, the **Heart® Patch** membrane acts as a matrix for fibroblast infiltration and as a substrate for the deposit of new collagen.

Heart® Patch is therefore **gradually degraded** and replaced by new, vital patient tissue.

Heart® Patch

completely **biocompatible**

hard-wearing and easy to **suture**

scaffold for the healing process


progressively replaced by **new vital tissue**



Heart® Patch is a natural membrane made from equine pericardium, treated using the sophisticated Zymo-Teck® deantigenation process


The natural matrix for the repair and strengthening of tendons and ligaments

The maintenance of the native pericardium structure, characterised by a weave of multi-directional, multi-layer collagen fibres makes **Heart® Patch** perfectly **adaptable** to the implant site and **high mechanical resistance** able to ensure easy fixing into place by suture.




Heart® Patch membrane observed under the scanning electronic microscope (SEM) has a multilayer, compact appearance, characterised by a close-knit weave of collagen fibres

Padua University, Biology Department, Electronic Microscopy Service



Recommended for:

- > Rotator cuff
- > Anterior shoulder capsule
- > Ulnar collateral ligament
- > Triceps brachii tendon
- > Flexors and extensors of the hand
- > Posterior hip capsule
- > Quadriceps tendon
- > Posterior tibial tendon
- > Peroneal tendon
- > Tendon of the knee
- > Achilles' tendon



Heart® Patch pericardium membrane

HRT-020	Pericardium membrane	50 x 50 x 0.2-0.4 mm
HRT-021	Pericardium membrane	60 x 80 x 0.2-0.4 mm
HRT-022	Pericardium membrane	∅ 50 x 0.2-0.4 mm

Application of **Heart® Patch** to **repair** and **strengthen** the posterior tibial tendon

Heart® Patch is cut to size and placed beneath the tendon to be strengthened



Heart® Patch is wound and sutured around the posterior tibial tendon to strengthen it and provide an excellent scaffold for the healing process



Photograph courtesy of Dr Alberto Siclari Trauma and Orthopaedic ward, Biella Hospital